

**Recombinant Bovine Fetuin A, Animal Component-Free**

<b>Cat. No. :</b>	B055C
<b>Alternative Names:</b>	Fetuin A; FETUA; Alpha-2-HS-Glycoprotein; Alpha-2-Z-Globulin; Ba-Alpha-2-Glycoprotein; Fetuin-A; AHSG; FETUA; α2-HS-glycoprotein
<b>Species:</b>	Bovine
<b>Accession No.:</b>	P12763
<b>Expression System:</b>	CHO
<b>Protein Sequence:</b>	aa 19-359
<b>Theoretical MW:</b>	36.35 kDa
<b>Theoretical pI:</b>	5.10
<b>Tag:</b>	Tag-Free.
<b>Formulation buffer:</b>	PBS, 5% Mannitol and 0.01% Tween 80, pH7.4.
<b>Appearance:</b>	Lyophilized Powder.
<b>Purity:</b>	≥95% as determined by SDS-PAGE.
<b>Bioactivity:</b>	Recombinant Bovine Fetuin A stimulates adhesion of B16-F1 cells, as measured in a cell proliferation assay. The ED <sub>50</sub> for this effect is ≤ 100 µg/mL.
<b>Endotoxin Level:</b>	≤0.01 EU/µg, as determined by the LAL assay.
<b>Application:</b>	Cell Culture; Activity Assays.

**Preparation & Storage**

<b>Reconstitution:</b>	Reconstitute with sterile double-distilled water (ddH <sub>2</sub> O). <div style="border: 1px solid orange; padding: 5px; margin-top: 5px;"> <p>⚠ Centrifuge the vial briefly before opening to ensure full recovery of the solution. Avoid vortexing and minimize vigorous pipetting to maintain protein stability.</p> <p>❄ Immediately aliquot the reconstituted protein solution and store under recommended conditions. Avoid repeated freeze-thaw cycles.</p> </div>
<b>Shipping:</b>	Shipped on dry ice. Short-term transit on cold packs (2-8°C) is acceptable.
<b>Storage:</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>● 12 months from date of receipt, -20 to -80°C as supplied.</li> <li>● 2-7 days at 2 to 8°C under sterile conditions after reconstitution.</li> <li>● 3-6 months at -20 to -80°C under sterile conditions after reconstitution.</li> </ul>

**Protein Description**

**Background:** Bovine Fetuin A (α2-HS-glycoprotein), encoded by AHSG, is a hepatocyte-derived 48-50 kDa serum glycoprotein (~ 1 mg/mL in bovine plasma). Structurally homologous to human fetuin A (~ 75% identity), it comprises two cystatin-like domains and a histidine-rich C-terminus. Functionally indispensable as:

- Calcification inhibitor: Forms soluble calciprotein particles (CPPs) with calcium/phosphate, preventing ectopic mineralization; the standard reagent for in vitro vascular calcification models.
- Cell culture cornerstone: Xeno-compatible additive in serum-free media (hybridoma, CHO, stem cells) to enhance attachment, reduce shear stress, and stabilize growth factors.
- Protease regulator: Cystatin-superfamily member with cysteine protease inhibitory activity.
- Physiological marker: Negative acute-phase protein in cattle (serum levels decline during inflammation).

Commercially purified bovine fetuin A is a gold-standard bioreagent due to high yield, stability, and cross-species functional conservation – bridging basic mineralization research, bioprocessing, and comparative physiology.

**References:**

1. Hermann BP, et al. Purification and characterization of bovine fetuin. J Biol Chem. 1984;259(15):9842-9848.
2. Müller D, et al. Molecular cloning and expression of bovine fetuin cDNA. Biochim Biophys Acta. 1990;1048(2):157-163.
3. Heiss A, et al. Structural basis of calcification inhibition by fetuin-A. J Biol Chem. 2005;280(14):13313-13318.
4. Jahnhen-Dechent W, Ketteler M. Fetuin-A (AHSG): a multifunctional protein. Nephrol Dial Transplant. 2005;20(10):2023-2027.
5. Freshney RI. Culture of Animal Cells: A Manual of Basic Technique and Specialized Applications. 7th ed. Wiley-Blackwell; 2016.
6. Speer MY, et al. Role of fetuin-A in vascular calcification. Pediatr Nephrol. 2010;25(11):2213-2219.

## Product Disclaimer

- This product is intended for research use only and is not for use in diagnostic, therapeutic, or clinical procedures.
- Due to the inherent nature of research reagents, any liability arising from product quality issues is limited solely to the replacement or refund of the product itself and excludes all indirect, incidental, or consequential damages.



获取前沿资讯、最新产品及活动

ANEWBiosystems 只为基于科学研究的用途提供服务，不为任何个人或者非科研性质的使用等其他用途提供服务

Tel: 86 + 132-0209-8186; 86 + 186-2016-1769 E-mail: [info@anewbiosystems.com](mailto:info@anewbiosystems.com)

A Foundation Built for Discovery. --- 探索之路，基石为筑

购买产品将视为接受 ANEWBiosystems 的销售条款及条件