

Recombinant Human Interleukin-1 beta (IL-1β), Animal Component-Free

Cat. No. :	H004E
Alternative Names:	IL1b; IL-1β; IL1β; IL-1B; IL1B; IL1F2; IL-1 beta; IL-1b; Interleukin-1 beta; Interleukin 1 beta; Catabolin
Species:	Human
Accession No.:	P01584
Expression System:	E. coli
Protein Sequence:	Ala117-Ser269
Theoretical MW:	17.38 kDa
Theoretical pI:	5.91
Tag:	Tag-Free.
Formulation buffer:	PBS, 5% Trehalose and 0.01% Tween 80, pH7.4.
Appearance:	Lyophilized Powder.
Purity:	≥95% as determined by SDS-PAGE.
Bioactivity:	The ED ₅₀ was ≤10 pg/mL, as measured by its ability to induce NF-κB reporter gene activity in MCF-7 cells.
Endotoxin Level:	≤0.1 EU/μg, as determined by the LAL assay.
Application:	Cell Culture; Activity Assays.

Preparation & Storage

Reconstitution:	<p>Reconstitute with sterile double-distilled water (ddH₂O).</p> <p>⚠ Centrifuge the vial briefly before opening to ensure full recovery of the solution. Avoid vortexing and minimize vigorous pipetting to maintain protein stability.</p> <p>❄ Immediately aliquot the reconstituted protein solution and store under recommended conditions. Avoid repeated freeze-thaw cycles.</p>
Shipping:	Shipped on dry ice. Short-term transit on cold packs (2-8°C) is acceptable.
Storage:	<p>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -80°C as supplied. ● 2-7 days at 2 to 8°C under sterile conditions after reconstitution. ● 3-6 months at -20 to -80°C under sterile conditions after reconstitution.

Protein Description

Background: Interleukin-1 beta (IL-1β) is a prototypical pro-inflammatory cytokine encoded by the IL1B gene in humans and is a central mediator of innate immunity, inflammation, and host defense. Unlike most secreted proteins, IL-1β lacks a conventional signal peptide and is synthesized as an inactive 31-kDa precursor (pro-IL-1β). Its biological activity requires proteolytic cleavage by caspase-1 – activated within the NLRP3 inflammasome in response to pathogen-associated or damage-associated molecular patterns (PAMPs/DAMPs)–to generate the mature, secreted 17-kDa form.

IL-1β exerts its effects by binding to the IL-1 receptor type I (IL-1RI), which recruits the IL-1 receptor accessory protein (IL-1RAcP) to form a signaling-competent complex. This triggers downstream activation of NF-κB, MAPK (e.g., p38, JNK), and AP-1 pathways, leading to the expression of inflammatory mediators such as IL-6, IL-8, TNF-α, COX-2, and adhesion molecules. Due to its potent bioactivity, IL-1β signaling is tightly regulated at multiple levels—including transcriptional control, inflammasome activation, and extracellular inhibition by the natural antagonist IL-1Ra (interleukin-1 receptor antagonist).

IL-1β plays critical roles in both physiological and pathological contexts:

- It is essential for host defense against infections but also drives chronic inflammation in diseases such as rheumatoid arthritis, gout, atherosclerosis, and type 2 diabetes.
- In Alzheimer’s disease, IL-1β is upregulated in activated microglia surrounding amyloid-β plaques and contributes to neuroinflammation and neuronal dysfunction.
- In cancer, IL-1β is overexpressed in multiple solid tumors (e.g., melanoma, lung, breast, colorectal, and head and neck cancers), where it promotes tumor growth, angiogenesis, metastasis, and immunosuppression. Elevated IL-1β correlates with poor prognosis and cancer-related cachexia.
- It also contributes to intervertebral disc degeneration and pulmonary fibrosis in experimental models.

Given its central role in inflammation and disease, IL-1β is a well-validated therapeutic target. Clinically approved agents such as anakinra (recombinant IL-1Ra), canakinumab (anti-IL-1β monoclonal antibody), and riloncept (IL-1 trap) are used to treat autoinflammatory syndromes and are being investigated in cardiovascular disease and oncology.

References:

1. Dinarello, C. A. (2011). Immunological and inflammatory functions of the interleukin-1 family. *Annual Review of Immunology*, 29, 707-742.
2. Latz, E., Xiao, T. S., & Stutz, A. (2013). Activation and regulation of the inflammasomes. *Nature Reviews Immunology*, 13(6), 397-411.
3. Rider, P., et al. (2011). IL-1 α and IL-1 β are endogenous ligands for the NLRP3 inflammasome. *European Journal of Immunology*, 41(12), 3471-3479.
4. Heneka, M. T., et al. (2018). Neuroinflammation in Alzheimer's disease. *The Lancet Neurology*, 14(4), 388-405. [https://doi.org/10.1016/S1474-4422\(15\)70016-5](https://doi.org/10.1016/S1474-4422(15)70016-5)
5. Voronov, E., et al. (2014). IL-1 in chronic inflammation: The link between obesity and cancer? *Cytokine & Growth Factor Reviews*, 25(6), 675-687.
6. Marchetti, C., et al. (2018). The IL-1 β pathway in inflammation-induced skeletal muscle wasting. *Journal of Cachexia, Sarcopenia and Muscle*, 9(4), 680-693.
7. Zhang, Y., et al. (2020). IL-1 β promotes intervertebral disc degeneration by inducing inflammatory responses in nucleus pulposus cells. *Frontiers in Immunology*, 11, 589654.

Product Disclaimer

- This product is intended for research use only and is not for use in diagnostic, therapeutic, or clinical procedures.
- Due to the inherent nature of research reagents, any liability arising from product quality issues is limited solely to the replacement or refund of the product itself and excludes all indirect, incidental, or consequential damages.



获取前沿资讯、最新产品及活动

ANEWBiosystems 只为基于科学研究的用途提供服务，不为任何个人或者非科研性质的使用等其他用途提供服务

Tel: 86 + 132-0209-8186; 86 + 186-2016-1769 E-mail: info@anewbiosystems.com

A Foundation Built for Discovery. --- 探索之路，基石为筑

购买产品将视为接受 ANEWBiosystems 的销售条款及条件