

**Recombinant Human Interleukin-3 (IL-3), Animal Component-Free**

<b>Cat. No. :</b>	H006C
<b>Alternative Names:</b>	IL3; IL-3; Interleukin-3; Interleukin 3; Hematopoietic growth factor; Mast cell growth factor; MCGF; Multipotential colony-stimulating factor; P-cell-stimulating factor
<b>Species:</b>	Human
<b>Accession No.:</b>	P08700
<b>Expression System:</b>	CHO
<b>Protein Sequence:</b>	Ala20-Phe152
<b>Theoretical MW:</b>	15.09 kDa
<b>Theoretical pI:</b>	7.05
<b>Tag:</b>	Tag-Free.
<b>Formulation buffer:</b>	PBS, 5% Mannitol and 0.01% Tween 80, pH7.4.
<b>Appearance:</b>	Lyophilized Powder.
<b>Purity:</b>	≥95% as determined by SDS-PAGE.
<b>Bioactivity:</b>	ED <sub>50</sub> ≤ 0.1 ng/mL, as determined in a TF-1 human erythroleukemia cell proliferation assay.
<b>Endotoxin Level:</b>	≤0.01 EU/μg, as determined by the LAL assay.
<b>Application:</b>	Cell Culture; Activity Assays.

**Preparation & Storage**

<b>Reconstitution:</b>	Reconstitute with sterile double-distilled water (ddH <sub>2</sub> O). <div style="border: 1px solid orange; padding: 5px; margin-top: 5px;"> <p>⚠ Centrifuge the vial briefly before opening to ensure full recovery of the solution. Avoid vortexing and minimize vigorous pipetting to maintain protein stability.</p> <p>❄ Immediately aliquot the reconstituted protein solution and store under recommended conditions. Avoid repeated freeze-thaw cycles.</p> </div>
<b>Shipping:</b>	Shipped on dry ice. Short-term transit on cold packs (2-8°C) is acceptable.
<b>Storage:</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>● 12 months from date of receipt, -20 to -80°C as supplied.</li> <li>● 2-7 days at 2 to 8°C under sterile conditions after reconstitution.</li> <li>● 3-6 months at -20 to -80°C under sterile conditions after reconstitution.</li> </ul>

**Protein Description**

**Background:** Interleukin-3 (IL-3), also known as multi-CSF (multilineage colony-stimulating factor), is a hematopoietic cytokine primarily secreted by activated CD4<sup>+</sup> T cells (especially Th2 subsets), mast cells, and NKT cells. It supports the survival, proliferation, and differentiation of multipotent hematopoietic stem and progenitor cells into myeloid lineages—including granulocytes, monocytes, erythrocytes, megakaryocytes, eosinophils, basophils, and mast cells.

IL-3 signals through a heterodimeric receptor composed of a specific α-subunit (CD123) and a shared β-subunit (CD131, common with IL-5 and GM-CSF receptors). Downstream signaling involves JAK2/STAT5, PI3K/Akt, and Ras/MAPK pathways, promoting cell growth and anti-apoptotic responses.

Beyond its role in steady-state and emergency hematopoiesis, IL-3 contributes to allergic inflammation (e.g., asthma) by enhancing basophil and mast cell activation. Aberrant IL-3 signaling—often via CD123 overexpression—is implicated in acute myeloid leukemia (AML), chronic myeloid leukemia (CML), and blastic plasmacytoid dendritic cell neoplasm (BPDCN), making CD123 a key therapeutic target.

Recombinant IL-3 is used in research for ex vivo expansion of hematopoietic progenitors, while clinical efforts focus on CD123-directed therapies (e.g., tagraxofusp, anti-CD123 CAR-T).

**References:**

1. Miyajima, A., et al. (1992). Interleukin-3 and its receptor. *Immunology Today*, 13(6), 204-208.
2. Sanderson, S., et al. (2021). IL-3: A cytokine at the interface of hematopoiesis and immunity. *Frontiers in Immunology*, 12, 705782.
3. Testa, U., et al. (2019). CD123: A novel therapeutic target in acute myeloid leukemia. *Cells*, 8(10), 1230.
4. Hermans, M. H., et al. (2003). Signaling by IL-3, IL-5, and GM-CSF: Common and unique elements. *Journal of Leukocyte Biology*, 73(6), 673-681.
5. Tagraxofusp (SL-401) FDA Approval Summary (2018). U.S. Food and Drug Administration.

## Product Disclaimer

- This product is intended for research use only and is not for use in diagnostic, therapeutic, or clinical procedures.
- Due to the inherent nature of research reagents, any liability arising from product quality issues is limited solely to the replacement or refund of the product itself and excludes all indirect, incidental, or consequential damages.



获取前沿资讯、最新产品及活动

ANEWBiosystems 只为基于科学研究的用途提供服务，不为任何个人或者非科研性质的使用等其他用途提供服务

Tel: 86 + 132-0209-8186; 86 + 186-2016-1769 E-mail: [info@anewbiosystems.com](mailto:info@anewbiosystems.com)

A Foundation Built for Discovery. --- 探索之路，基石为筑

购买产品将视为接受 ANEWBiosystems 的销售条款及条件